

Superior Technology

Blockparty®

Innovation



Service

Value

Imagine being able to house mice in up to 750 square inches of open caging! Easy cage-to-cage connectivity on Optimice® rack up to 10 in all per rack level.

Researchers can now study mice in an environment that encourages natural behavior:

- Exploration
- Preference testing
- Social grouping and hierarchy identification
- Territory and resource competition
- Exercise and fitness
- Breeding
- Nesting

Available with amber or red tube assembly

Improve the
Productivity
of your Lab
Today!

All of our products and custom accessories are listed and described in our catalog at: www.animalcaresystems.com